THUASNE®

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en LONG-STRETCH ELASTIC COMPRESSION BANDAGE

Description/Destination

The product comes in two versions (16 and 17), offering different levels of pressure depending on the coverage used ①.

Pressure 1: 20 to 45 cN/cm - **Pressure 2:** 46 to 100 cN/cm - **Pressure 3:** 101 to 160 cN/cm - **Pressure 4** (for information only): > 160 cN/cm.

Each of the bandages is available in different versions (graduated or non-graduated, with and without stirrup) and different widths, depending on need. Biflex is intended for use on the upper and lower limbs of adult or pediatric patients who meet one of the defined medical indications.

Composition: polyamide - viscose - polyester - elastane.

Properties/Mode of action

Biflex bandages are long-stretch two-way bandages:

longitudinal stretch > 100%, springback > 90%, transverse elasticity > 40%.

Indications

Alternative to compression stockings (15-36 mmHg) for varicose veins. Chronic edema. Post-traumatic edema. After sclerotherapy or surgery of varicose veins as an alternative to compression stockings. Alternative to compression stockings for trophic disorders (dermatitis/lipodermatosclerosis). Venous leg ulcers (multitype bandaging, elastic bandage for patients with reduced mobility). Symptomatic treatment of deep or superficial acute venous thrombosis. Alternative to compression or antithrombotic stockings. Prevention of post-thrombotic syndrome: alternative to compression stockings for at-risk patients. Treatment of lymphedema.

Contraindications

Do not use in case of major dermatological conditions on an affected limb. Do not use in the event of known allergy to any of the components. Do not use in case of severe peripheral neuropathy of the affected limb. Do not use in the event of septic thrombosis. Do not use in the event of phlegmasia cerulea dolens (painful blue phlebitis with arterial compression) of an affected limb.

Specific to the lower limbs:

Do not use in case of lower limb peripheral arterial disease (PAD) with ankle-brachial index (ABI) < 0.6. Do not use in the event of decompensated heart failure. Do not use in the event of advanced diabetic microangiopathy (for a compression > 30 mmHg). Do not use in case of extra-anatomical bypass of the affected limb.

Specific to the upper limbs:

Do not use in case of brachial plexus condition.
Do not use in case of vasculitis of the extremities.

Precautions

Strictly comply with your healthcare professional's prescription and recommendations for use. The risk/benefit ratio and the adequate level of pressure must be regularly re-assessed by a healthcare professional in the event of:

- PAD with ABPI between 0.6 and 0.9
- Advanced peripheral neuropathy
- Weeping or eczematous dermatitis.

In the event of discomfort, significant restriction, pain, variation in limb

volume, deterioration of skin condition, infection, unusual sensations, change in the colour of the extremities, or changes in the performance of the product, remove the product and seek the advice of a healthcare professional. Do not apply the product directly to an open wound without a dressing. For hygiene, security and performance reasons, do not re-use the product for another patient. Do not use products on the skin (creams, ointments, oils, gels, patches...) before applying the device, it could damage the product. This product can be combined with other bandages for use during the day and at night. Consult a healthcare professional and follow the instructions if you need to remove the bandage at night. The choice of the width is at the discretion of the healthcare professional, depending on the limb to be treated. In the event of a venous ulcer, a Doppler examination is necessary before using the device in order to ensure there is no associated severe arterial disease. For some conditions (or situations) such as venous thrombosis, the device should be used in combination with an anticoagulant treatment; consult a healthcare professional for advice.

Undesirable side-effects

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity or cutaneous dryness. Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

Check that the product is undamaged before use.

Examine the shape of the leg; if there are any deformed or bony areas, use protective foam bandages or Varico pads to make the shape of the limb uniform.

For a lower limb:

1. Position the end of the bandage flat against the base of the toes.
2. Wind round a first time going under the foot.
3. The second wind should go towards the ankle. Cover the heel.
4. Go under the foot and back to the ankle.
5. Return to the top of the foot and go round the ankle.

6. Go up the leg, spacing the turns regularly, in accordance with your prescription.
7. Close the bandage with adhesive tape or a special fastener.

For an upper limb:

1. Position the end of the bandage flat against the base of the fingers.
2. Make one or two loops around the hand to hold the bandage in place, if necessary.
3. Go up the arm, spacing the turns regularly, in accordance with your prescription.
4. Fix the end of the bandage in place with sticking plaster or a dedicated fastening system.

If using graduated bandages:

The markings help control the level of pressure when applying the bandage. The bandage is correctly stretched when each rectangle is transformed into a square ❶. The markings also help to guide the coverage:

- Half the square covered = pressure 1 (BIFLEX 16+)/pressure 3 (BIFLEX 17+) ❷.
- Full square covered = pressure 2 (BIFLEX 16+)/pressure 4 (BIFLEX 17+) ❸.

NB: the level of coverage is to be determined by the healthcare professional when using non-graduated bandages.

Biflex + Pratic or stirrup:

This bandage has a loop at one end to attach it to the foot, as well as a transfer indicating which way the bandage is to be wound.

NB: when the bandage is removed, it needs to be unwound from the end without the loop.

Care

Wash before every use. Machine washable at 40°C (delicate programme). Do not use detergents, fabric softeners or aggressive products (products containing chlorine, etc.). Squeeze out excess water. Dry away from direct sources of heat (radiator, sun, etc.). Do not bleach. Do not tumble-dry. Do not iron. Do not dry clean. Can be autoclaved.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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©Thuasne - 2047501 (2023-04)
Thuasne UK Ltd
Thuasne SAS – SIREN/RCS Nanterre 542 091 186 – capital 1 950 000 euros