THUASNE

BIFLEX 16

BIFLEX 17

BIFLEX 16+

BIFLEX 17+

BIFLEX 16+ PRATIC


BIFLEX 17+ PRATIC

MD

fr	Dispositif médical
en	Medical device
de	Medizinprodukt
nl	Medisch hulpmiddel
it	Dispositivo medico
es	Producto sanitario
pt	Dispositivo médico
da	Medicinsk udstyr
fi	Lääkinnällinen laite
sv	Medicinteknisk produkt
el	Ιατροτεχνολογικό προϊόν
pl	Wyrób medyczny
lv	Medicīniska ierīce
sk	Zdravotnícka pomôcka
hu	Orvostechnikai eszköz
bg	Медицинско изделие
ro	Dispozitiv medical
ru	Медицинское изделие
hr	Medicinski proizvod
zh	醫療器械
ja	医療機器
ko	의료 기기
ar	جهاز طبي


fr	Un seul patient - à usage multiple
en	Single Patient - multiple use
de	Einzelner Patient - mehrfach anwendbar
nl	Eén patiënt - meervoudig gebruik
it	Singolo paziente - uso multiplo
es	Un solo paciente - uso múltiple
pt	Paciente único - várias utilizações
da	Enkelt patient - flergangsbrug
fi	Potilaskohtainen - voidaan käyttää useita kertoja
sv	En patient - flera användningar
el	Πολλαπλή χρήση - σε έναν μόνο ασθενή
pl	Wielokrotne użycie u jednego pacjenta
lv	Viens pacients - vairākkārtēja lietošana
sk	Jeden pacient - viacnásobné použitie
hu	Egyetlen beteg esetében többször újrahasználható
bg	Един пациент - многократна употреба
ro	Un singur pacient - utilizare multiplă
ru	Многократное использование для одного пациента
hr	Jedan pacijent - višestruka uporaba
zh	一位患者 - 多次使用
ja	一患者用 - 複数回使用
ko	한 환자에 한 해 - 여러번 사용
ar	لمريض واحد - متكرر الاستعمال

BIFLEX 16**BIFLEX 17**

fr	Comment différencier les bandes BIFLEX: 1 fil bleu = BIFLEX 16 - 2 fils bleus = BIFLEX 17
en	How to identify BIFLEX bandages: 1 blue thread = BIFLEX 16 - 2 blue threads = BIFLEX 17
de	Unterscheidungsmerkmale der verschiedenen BIFLEX-Binden: 1 blauer Faden = BIFLEX 16 - 2 blaue Fäden = BIFLEX 17
nl	Hoe herkent u de verschillende pleisters van BIFLEX: 1 blauwe draad = BIFLEX 16 - 2 blauwe draden = BIFLEX 17
it	Come distinguere le bende BIFLEX: 1 filo blu = BIFLEX 16 - 2 fili blu = BIFLEX 17
es	Cómo diferenciar las vendas BIFLEX: 1 hilo azul = BIFLEX 16 - 2 hilos azules = BIFLEX 17
pt	Como distinguir as ligaduras BIFLEX: 1 fio azul = BIFLEX 16 - 2 fios azuis = BIFLEX 17
da	Sådan kender man forskel på BIFLEX-bindene: 1 blå tråd = BIFLEX 16 - 2 blå tråde = BIFLEX 17
fi	Kuinka erottaa BIFLEX-siteet toisistaan: 1 sininen lanka = BIFLEX 16 - 2 sinistä lankaa = BIFLEX 17
sv	Hur du ser skillnad på BIFLEX-bandage: 1 blå tråd = BIFLEX 16 - 2 blå trådar = BIFLEX 17

el	Πώς να ξεχωρίσετε τους επιδέσμους BIFLEX: 1 μπλε νήμα = BIFLEX 16 - 2 μπλε νήματα = BIFLEX 17
pl	Jak odróżnić rodzaje bandaży BIFLEX: 1 pojedynczy niebieski ścieg = BIFLEX 16 - 2 podwójne niebieskie ściegi = BIFLEX 17
lv	Kā atšķirt saites BIFLEX: 1 zils diegs = BIFLEX 16 - 2 zili diegi = BIFLEX 17
sk	Ako rozlíšiť obvazy BIFLEX: 1 modré vlákno = BIFLEX 16 - 2 modré vlákna = BIFLEX 17
hu	Hogyan különböztethetők meg egymástól a BIFLEX fáslik: 1 kék csík = BIFLEX 16 - 2 kék csík = BIFLEX 17
bg	Как да различаваме бинтовете BIFLEX: 1 синя нишка = BIFLEX 16 - 2 сини нишки = BIFLEX 17
ro	Cum diferentțiem benzile BIFLEX: 1 fir albastru = BIFLEX 16 - 2 fire albastre = BIFLEX 17
ru	Как различаются бинты BIFLEX: 1 голубая полоска – BIFLEX 16 - 2 голубые полоски – BIFLEX 17
hr	Kako razlikovati BIFLEX zavoje: 1 plava linija = BIFLEX 16 - 2 plave linije = BIFLEX 17
zh	如何区别BIFLEX医用绷带: 1条蓝线 = BIFLEX 16 - 2条蓝线 = BIFLEX 17
ja	BIFLEX (バイフレックス) 包帯の識別方法 青筋 1本 = BIFLEX 16 - 青筋 2本 = BIFLEX 17
ko	BIFLEX 붕대의 종류를 구별하는 방법: 1개의 파란색 선 = BIFLEX 16 - 2개의 파란색 선 = BIFLEX 17
ar	كيفية التفرقة بين أربطة BIFLEX: BIFLEX 17 = خيط واحد أزرق - BIFLEX 16 = خيطان أزرقان

1

2

3

fr

Description

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en

Description

Long-stretch elastic compression bandage

Composition:

polyamide, viscose, polyester, elastane

Properties/Mechanism of action

BIFLEX bandages are long-stretch bandages indicated for the treatment of venous and lymphatic conditions as well as trauma. The bandages are made of a two-way stretch fabric: longitudinal stretch > 100%, elastic return > 90%, transverse stretch > 40%.

Use:

BIFLEX is intended for use on the arms or legs and is suitable for adults and children with one of the listed medical indications.

Indications

Alternative to compression stockings (15-36 mmHg) for varicose veins

Chronic edema

Post-traumatic edema

After sclerotherapy or surgery of varicose veins as an alternative to compression stockings

Alternative to compression stockings for trophic disorders (dermatitis/lipodermatosclerosis)

Venous leg ulcers (multitype bandaging, elastic bandage for patients with reduced mobility)

Symptomatic treatment of acute superficial and deep venous thrombosis

Alternative to compression or antithrombotic stockings

Post-thrombotic syndrome prevention: alternative to compression stockings for high-risk patients

Lymphedema treatment

Fitting instructions

Check the stocking is intact before use.

1. Place one end of the bandage flat against the base of the toes.
2. Wind round once, passing the bandage under the foot.
3. Wind round a second time, up towards the ankle. Cover the heel.
4. Pass the bandage back under the foot and return to the ankle.

5. Return to the top of the foot and pass the bandage round the ankle.
6. Wind the bandage up round the calf, evenly spacing the rounds according to the instructions.

7. Secure the bandage using a sticking plaster or the fastener provided. The bandage may be used together with a foam bandage or Varico foam pads for deformed or bony areas.

Graduated bandages:

The markings help control the level of compression when applying the bandage. The correct level of stretch is achieved when the rectangles printed on the bandage are stretched into squares (diagram 1). The markings also act as a guide for the level of coverage:

- Squares half covered with each round = Level 1 (BIFLEX 16+) / Level 3 (BIFLEX 17+) (diagram 2)
- Squares fully covered with each round = Level 2 (BIFLEX 16+) / Level 4 (BIFLEX 17+) (diagram 3)

BIFLEX + PRATIC bandages:

The bandage has a stirrup at one end to secure it around the foot, with markings showing the direction in which to wrap the bandage.

NB: After removing the bandage, roll it up starting at the end without the stirrup.

Contraindications

PAD (peripheral arterial disease) with ABI < 0.6

Skin wounds, abscess, infected skin conditions

Extra-anatomic bypass

Uncontrolled heart failure

Phlegmasia cerulea dolens (painful blue oedema with arterial compression)

Septic thrombosis

Hypersensitivity to any of the components

Severe peripheral neuropathy

Precautions

Follow the advice of the professional who prescribed or supplied the product.

Examine the shape of the leg; if there are any deformed or bony areas, use protective foam or Varico pads to even out the surface of the limb.

In case of discomfort, significant impediment, pain, change in the size of the leg or unusual sensations, remove the device and consult a health

care professional.

Follow the instructions from your doctor or health care professional for how to fit the bandage and how long to wear it.

This device may be combined with other bandages for day-time and night-time use. If you need to remove the bandage at night, consult your doctor and follow the instructions given.

Store at room temperature, preferably in the original packaging.

Do not apply directly to an open wound.

Fully cover any infected or open skin wounds before using the bandage.

For some conditions (or situations) like venous thrombosis, bandages should be used in combination with anticoagulants; please refer to the instructions from your health care professional.

Any serious incidents relating to the device must be reported to the manufacturer and competent authority in the Member State in which the user and/or the patient is based.

Adverse effects

This device may cause a skin reaction such as pressure marks, redness, itching, burning or blisters.

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Polyamid, Viskose, Polyester, Elasthan

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©Thuasne - 2019501 (2020-02)