


THUASNE®

# LombaStab

fr	Ceinture de soutien lombaire - "Serrage facile et précis"*	5
en	Lumbar support belt - "Easy, precise tightening"*	6
de	LWS-Stabilisierungsorthese - „einfaches und präzises Justieren“*	7
nl	Lumbale rugbrace - "Eenvoudige en nauwkeurige aanspanning"*	8
it	Cintura di sostegno lombare - "Serraggio facile e preciso"*	9
es	Faja de sujeción lumbar - "Fijación fácil y precisa"*	11
pt	Cinta de suporte lumbar - "Aperto fácil e preciso"*	12
da	Rygstøttebælte - "Nem og præcis stramning"*	13
fi	Lannerangan tukivyo - "Helpostija täsmällisesti kiristettävä"*	14
sv	Stödbälte för ryggen - "enkel och precis åtdragnings"*	15
el	Zώνη οσφυϊκής υποστήριξης - «Εύκολη και ακριβής σύσφιξη»*	17
cs	Bederní páš - „Jednoduché a přesné stažení“*	18
pl	Orteza lędźwiowo-krzyżowa - „Łatwe i precyzyjne mocowanie“*	19
lv	Muguras lejasdājas josta - „Viegla un precīza savilkšana“*	20
lt	Palaikomasis juosmens diržas - „Lengvai ir tiksliai suveržiamas“*	21
et	Nimme piirkonna tugivõö - „Lihtne ja täpne pingutus“*	23
sl	Oporni križni pas - „Enostavno in natančno vpetje“*	24
sk	Podporný driekový páš - „jednoduché a presné dotiahnutie“*	25
hu	Ágyék fűző - „Egyeszerű és pontos rögzítés“*	26
bg	Лумбален колан - „лесно и точно затягане“*	27
ro	Centură de susținere lombară - „Strângere ușoară și precisă“*	28
ru	Опорный поясничный корсет - «Легко и просто затянуть»*.	30
hr	Lumbalni potporni pojaz - „Jednostavno i precizno stezanje“*	31
zh	易于调节的加强型腰部稳定护带*	32
ar	حزام لدعم أسفل الظهر - "شد سهل ودقيق".....	34

	1	60 - 80 cm
	2	75 - 95 cm
	3	90 - 110 cm
	4	105 - 125 cm
	5	120 - 140 cm

<b>fr</b>	Tour de taille	Hauteurs : 21 cm - 26 cm
<b>en</b>	Waist circumference	Heights: 21 cm - 26 cm
<b>de</b>	Taillenumfang	Höhen: 21 cm - 26 cm
<b>nl</b>	Tailleomvang	Hoogtes: 21 cm - 26 cm
<b>it</b>	Circonferenza vita	Altezze: 21 cm - 26 cm
<b>es</b>	Contorno de cintura	Alturas: 21 cm - 26 cm
<b>pt</b>	Perímetro de cintura	Alturas: 21 cm - 26 cm
<b>da</b>	Taljemål	Højde: 21 cm - 26 cm
<b>fi</b>	Vyötärönympärys	Korkeudet: 21 cm - 26 cm
<b>sv</b>	Midjemått	Höjd: 21 cm - 26 cm
<b>el</b>	Περίμετρος μέσης	Υψη: 21 εκ. - 26 εκ.
<b>cs</b>	Obvod pasu	Výšky: 21 cm - 26 cm
<b>pl</b>	Obwód w talii	Wysokości: 21 cm - 26 cm
<b>lv</b>	Vidukļa apkārtmērs	Augstumi: 21 cm - 26 cm
<b>lt</b>	Juosmens apimtis	Aukštis: 21 cm - 26 cm
<b>et</b>	Võõrumbermõõt	Kõrgused: 21 cm - 26 cm
<b>sl</b>	Obseg pasu	Višina: 21 cm - 26 cm
<b>sk</b>	Obvod pása	Výšky: 21 cm - 26 cm
<b>hu</b>	Derékbőség	Magasság: 21 cm - 26 cm
<b>bg</b>	Обиколка на талията	Височини: 21 см - 26 см
<b>ro</b>	Circumferința taliei	Înălțimi: 21 cm - 26 cm
<b>ru</b>	Обхват талии	Высота: 21 см - 26 см
<b>hr</b>	Opseg struka	Visine: 21 cm - 26 cm
<b>zh</b>	腰围	高度 : 21 cm -26 cm
<b>ar</b>	محيط الخصر	الارتفاعات: 21 سم - 26 سم

<b>fr</b>	Compression	Maintien	Éducation posturale
<b>en</b>	Compression	Support	Postural education
<b>de</b>	Kompression	Halt	Erlernen der Körperhaltung
<b>nl</b>	Compressie	Fixatie	Voorlichting over houding
<b>it</b>	Compressione	Sostegno	Educazione posturale
<b>es</b>	Compresión	Sujeción	Educación postural
<b>pt</b>	Compressão	Manutenção	Educação postural
<b>da</b>	Kompression	Støtte	Postural genoptræning
<b>fi</b>	Kompressio	Tuki	Asentokuntoutus
<b>sv</b>	Kompression	Stöd	Postural träning
<b>el</b>	Συμπίεση	Συγκράτηση	Σωστή στάση σώματος
<b>cs</b>	Komprese	Stabilizace	Správné držení těla
<b>pl</b>	Kompresja	Podtrzymywanie	Korekcja postawy
<b>lv</b>	Kompresija	Atbalsts	Stājas korekcija
<b>lt</b>	Kompresija	Palaikymas	Laikysenos taisymas
<b>et</b>	Kompressioon	Hooldus	Rühitreening
<b>sl</b>	Kompresija	Orijem	Učenje telesne drže
<b>sk</b>	Kompresia	Podpora	Cvičenie posturálnych svalov
<b>hu</b>	Szorítás	Támasz	Tartásjavító oktatás
<b>bg</b>	Компресия	Поддръжка	Приучаване към правилна стойка
<b>ro</b>	Compresie	Susținere	Educație posturală
<b>ru</b>	Компрессия	Стабилизация	Исправление осанки
<b>hr</b>	Kompresija	Potpora	Učenje o držanju
<b>zh</b>	压紧	稳定	姿势联系
<b>ar</b>	الضغط	الدعم	تدريب وضعی

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Conserver cette notice.

en

## LUMBAR SUPPORT BELT - "EASY, PRECISE TIGHTENING"\*\*

### Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table. This device can be used to support rehabilitation, particularly when resuming activities.

### Composition

Elastic fabric. ①

Handles. ②

Back panel with Coolmax lining. ③

Abdominal panels with finger-loops. ④

Laces.

Textile components: polyamide - polyester - elastane - viscose - polyurethane - polyethylene.

Non-textile components: brushed tempered steel - polypropylene - SEBS - polyamide - polyurethane - polyoxymethylene.

### Properties/Mode of action

The device relieves lower back pain by supporting the spine through an increase in intra-abdominal pressure. This reduces the load on the intervertebral discs.

This device uses back pads, implementing a pulley-like system to distribute pressure.

### Indications

Non-specific low back pain (acute, subacute and chronic phases).

Spinal stenosis.

## Contraindications

Do not use in the event of known allergy to any of the components.

Do not use on pregnant women.

Do not use in the event of a hiatus hernia.

Do not use in the event of bone cancer with spinal metastases.

Do not use in the event of circulatory, pulmonary or cardiovascular problems in patients for whom an increase in arterial blood pressure is not recommended.

## Precautions

Verify the product's integrity before every use.

Wear the device above thin clothing.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

Strictly comply with your healthcare professional's prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

Do not wear the device while sleeping.

In the event of discomfort, significant hindrance, pain, abnormal sensations, remove the device and consult a healthcare professional.

The use of a medical device by a child should be done under the supervision of an adult or a healthcare professional.

Measure the waist circumference of the patient at the mid-distance between the iliac crests and floating ribs.

It is recommended that a healthcare professional supervises the first application.

For hygiene, security and performance reasons, do not re-use the product for another patient.

Do not apply the product in direct contact with skin (risk of pinching).

For safety reasons, do not wear the belt close to machines (risk of becoming caught).

Do not wear the device in a medical imaging machine.

## Undesirable side-effect

This device can cause skin reactions (redness, itching, irritations...).

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

## Instructions for use/Application

It is recommended that a healthcare professional adjust the stays to the shape of the patient's back before the belt is worn for the first time.

Open the belt and lay it out completely flat so that you can see the handles. ②

Check that the two tightening handles ② are fully retracted.

If necessary, unhook them and separate the two back parts by pulling the ends of the belt, then fasten the handles again.

Position the belt around the waist, using the finger loops located on the front panels. ④

The tightening back pads must be positioned outwards (i.e. : not against the skin).

The lacing section must be centered over the spine: the pads should be positioned on either side of the spine. ⑤

The bottom of the belt must stop at the top of the crease between the buttocks. ⑥

The Thuisnies® logo indicates which way up the belt should be put on.

Fasten the belt: right side over left side, without pulling too tight or applying any particular pressure. ⑦ ⑧

Adjust the tightness using the two handles: pull them forwards then fold them flat against the front of the belt, where they attach like a self-fastener. ⑨ ⑩

These two handles can be used to adjust the level of compression as required throughout the day.

#### Care

Close the self-fastening tabs before washing. Machine washable at 30°C (delicate programme). Do not machine wash the device more than 10 times. Do not dry clean. Do not use detergents, fabric softeners or aggressive products (products containing chlorine, etc.). Do not tumble-dry. Do not iron. Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat.

#### Storage

Store at room temperature, preferably in the original packaging.

#### Disposal

Dispose of in accordance with local regulations.

\*Study conducted in-house on a panel of 13 people, July 2019

This medical device is a regulated health product with CE marking.

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©Thuasne - 2054001 (2024-04)
**THUASNE UK Ltd**
Thuasne SAS – SIREN/RCS Nanterre 542 091 186