
THUASNE®

en	Sport open strapping knee brace	4
fr	Genouillère strapping ouverte sport	5
de	Offenene Strapping Sport-Kniebandage	7
nl	Open kneibrace met elastische banden sport	9
it	Ginocchiera strapping aperta per uso sportivo	11
es	Rodillera strapping abierta para uso deportivo	13
pt	Joelheira de "strapping" aberta para uso desportivo	15
da	Åben knæbandage sport	17
fi	Nauhallinen avoin polvituki urheilukäytöön	18
sv	Knästöd med öppet elastiskt band sport	20
el	Αθλητική επιγονατίδα strapping ανοικτή	22
cs	Otevřená kolenní ortéza strapovací sport	24
pl	Bandażowa ortezka kolanowa otwarta sportowa	26
lv	Atvērta sporta elastīga ceļa ortoze	28
lt	Apjuosiamas atviras kelio įtvaras sportui	29
et	Spordi avatud põlvepealsega põlvekaitse	31
sl	Odprta opornica iz trakov za koleno šport	33
sk	Športová otvorená kolenná ortéza s popruhmi	34
hu	Nyitott pántos térdmerevítő sportoláshoz	36
bg	Отворена наколенка с каишка за спорт	38
ro	Genunchieră deschisă cu benzi ajustabile sport	40
hr	Otvoreni sport steznik za koljeno s trakama	42

SPORT

Instructions for use - Instructions d'utilisation - Gebrauchsanweisung - Gebruiksaanwijzingen - Istruzioni per l'uso - Instrucciones de uso - Instruções de utilização - Anvisninger vedrørende brug - Käyttöohjeet - Råd vid användning - Οδηγίες χρήσης - Pokyny k používání - Szczegółowe użycowania - Lietošanas instrukcija - Naudojimo instrukcijos - Kasutusjuhised - Navodila za uporabo - Pokyny na používanie - Használati utasítás - Инструкции за употреба - Instrucțiuni de utilizare - Upute za uporabu

	cm
S	32 - 36
M	37 - 40
L	41 - 44
XL	45 - 48

en	Knee circumference in slight flexion
fr	Circonférence de genou en légère flexion
de	Knieumfang bei leichter Beugung
nl	Omvang van de licht geplooide knie
it	Circonferenza del ginocchio in leggera flessione
es	Circunferencia de la rodilla con ligera flexión
pt	Circunferência do joelho em ligeira flexão
da	Knæets omkreds i let bøjning
fi	Polven ympärys hieman koukitettuna
sv	Omkrets runt knäet i lätt böjning
el	Περιφέρεια γόνατος σε ελαφρά κάμψη
cs	Obvod kolene v mírném ohýbu
pl	Obwód kolana przy lekkim zgięciu
lv	Celgala apkārtmērs nelielā saliekumā
lt	Šiek tiek sulenkto kelio apimtis
et	Kergelt pinguldatult põlve ümber
sl	Obseg kolena pri rahlem upogibu
sk	Obvod kolena s miernym ohybom
hu	Térd kerülete enyhe hajlításban
bg	Обиколка на коляното при лека флексия
ro	Circumferința genunchiului în flexiune ușoară
hr	Obuhvaća koljeno u blagoj fleksiji

SPORT OPEN STRAPPING KNEE BRACE

Ligament knee brace with functional strapping

Support level 3

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table. Bilateral product.

Available in 4 sizes.

The device is composed of:

- flexible lateral reinforcement,
- a double functional strapping ensuring good lateral support,
- a synthetic rubber (SBR) tissue.

This device can be used when resuming sports activities.

Composition

Textile components: polyamide - synthetic rubber (SBR) - polyester - thermoplastic polyurethane - elastane - ethylene vinyl acetate - polyethylene - natural latex.

Rigid components: polyamide - zinc alloy.

Properties/Mode of action

Stabilisation of the knee ligaments provided by the double strapping system.

Anti-slip system to prevent the knee brace slipping.

Opening to make fitting easier.

Limitation of the heat loss by the synthetic rubber (SBR).

Indications

Joint instability/laxity.

Knee pain and/or swelling management.

Contraindications

Do not use in the event of known allergy to any of the components.

Do not apply the product in direct contact with broken skin.

Do not use in the event of major venous thromboembolic history without thromboprophylaxis.

Precautions

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

For hygiene and performance reasons, do not re-use the product for another patient.

Consult a healthcare professional in case of uncertain diagnosis.

It is recommended to adequately tighten the device to achieve support without restricting blood circulation.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

In case of history of venous or lymphatic disorders, consult a healthcare professional.

In case of doubt before any sports activity, consult a healthcare professional.

Contains natural latex.

In the event of any modification in the product's performance, remove it and consult a healthcare professional.

Do not wear the device in a medical imaging machine.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters, etc.) or wounds of various degrees of severity.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

It is recommended to wear the device directly on the skin, unless contraindicated.

To avoid damage to the knit fabric, please make sure that the self-gripping pads are always fastened to the designated self-gripping surfaces when putting it on and taking it off.

Close the tissue fabrics sheath onto the joint as slightly bent, starting from the low portion. Adjust the kneecap adjustment first ① and then the fastening tab ②. Repeat the operations for the upper portion ③.

Once in place, the knee brace clamps and holds the kneecap.

Tighten the lower cross-over straps and wrap them around the rear of the knee to fasten them ④ ⑤.

Then tighten the upper straps and fasten them at rear of the knee support ⑥ ⑦.

Make sure the straps do not go on the kneecap. Release it if necessary. For comfort during wear, make sure the straps are not creased.

If you feel the need to tighten the straps during an activity, make sure you tighten them following the steps indicated in the fitting instructions.

Care/Maintenance

Close the self-fastening tabs before washing. Machine washable at 30 °C. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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©Thuasne - 2049701 (2023-10)
Thuasne UK Ltd
Thuasne SAS – SIREN/RCS Nanterre 542 091 186